

Limited Warranty

Wellness Arch is warranted against defects in wood body workmanship for 5 years after purchase.

Electronic controller is warranted for 2 years after purchase.

The shipping costs for replacement parts are waived during the first year after purchase. Additionally, all the shipping costs for returning parts or the device itself must be paid by the purchaser from one year after purchase.

The warranty does not apply to damages from accident, misuse, or alteration; nor does it apply to equipment being improperly operated or wired or maintained in violation of the unit's operation instructions.

All warranty claims must be accompanied by proof of the date of the purchase.

Manufacturer
Arch Tech Inc.
www.wellness-arch.com

User Manual



Wellness Arch

Handcrafted radiant heat therapeutic device

Arch Tech Inc.

Wellness Arch

Safety Guide for Electric Power

- This unit is designed to operate on AC current only.
- Never touch the power plug with wet hands.
- To prevent damage to the power cord, do not bend excessively or allow weighted objects to be placed on it.
- Never use a damaged power cord or outlet.
- Simultaneous use of other appliances on the same circuit may blow fuses or cause the unit to malfunction.
- Dusty power plug terminals should be wiped with dry cloth before use.
- If any smoking or burning smell should emit from the unit, please stop using it immediately and remove the power cord from the outlet. Contact the manufacturer.

Safety Guide for Maintenance

- Please do not place any flammable material near the unit.
- For cleaning the unit, please use dry cloth.
- Do not spray water on unit. Do not use benzene, thinner or bleach for cleaning the unit.

Product Handling **IMPORTANT**

- Please do not touch the heated area while in use. Injury may occur. Always grab edges where there is no heat for moving.
- Always use the unit at comfortable temperatures. Do not use at excessively hot temperatures. Hotter temperature does not provide better therapeutic results!

How to Operate

1. Turn on power by pressing power button.
2. Set the timer by pressing Timer Button if necessary. "0" indicates continuous operation.
3. Set the temperature using up and down buttons.

How to Use

Use the Wellness Arch 10 to 30 minutes per session for local therapy. You may increase the duration of use as you experience the device. Always start from lower temperatures. Below are examples of use:

