

Limited Warranty

Wellness Arch mini is warranted against defects in wood body workmanship for 5 years after purchase.

Electronic controller is warranted for 2 years after purchase.

The shipping costs for replacement parts are waived during the first year after purchase. Additionally, all the shipping costs for returning parts or the device itself must be paid by the purchaser from one year after purchase.

The warranty does not apply to damages from accident, misuse, or alteration; nor does it apply to equipment being improperly operated or wired or maintained in violation of the unit's operation instructions.

All warranty claims must be accompanied by the proof of the date of the purchase.

Manufacturer
Arch Tech Inc.
www.wellness-arch.com

User Manual



Wellness Arch Mini

Handcrafted radiant heat therapeutic device

Arch Tech Inc.

Wellness Arch

How to Operate



To turn on: Press any button

To turn off: Press the button in the middle for a few seconds

To increase temperature: Press the left button until the desired temperature is reached. Stay pressing the button to reach the top temperature. The indication light will get brighter.

To reduce temperature: Press the button in the middle until the desired temperature is reached. The indication light will be dimmed.

Right button changes the mode of the operation. You may use it for fun to enjoy the various indication light performances. To go back to the standard mode, press the left button.

How to Use

Use the Wellness Arch mini 10 to 30 minutes per session for local therapy. You may increase the duration of use as you experience the device. Always start from the lower temperature. Below are examples of use:

- Soothe tired eyes
- Promote blood circulation around scalp
- Relieve muscle and joint pain anywhere
- Use on upset stomach for stomach pain, diarrhea, and menstrual cramps
- Use as a bed warmer
- Temporarily reduce tooth pain
- Soothe cold knees
- improve skin condition by rubbing face or other areas

Safety Guide for Electric Power

- Never touch the power plug with wet hands.
- To prevent damage to the power cord, do not bend excessively or allow weighted objects to be placed on it.
- Never use a damaged power cord or outlet.
- Dusty power plug terminals should be wiped with dry cloth before use.
- If any smoking or burning smell should emit from the unit, please stop using it immediately and remove the power cord from the outlet. Contact the manufacturer.

Safety Guide for Maintenance

- Please do not place any flammable material near the unit.
- For cleaning the unit, please use dry cloth.
- Do not spray water on unit. Do not use benzene, thinner or bleach for cleaning the unit.

Product Handling **IMPORTANT**

- Always use the unit at a comfortable temperature. Do not use at excessively hot temperatures to avoid burning and/or injury. Long term use may result in burning and/or injury even at lower temperatures and requires extra caution. Hotter temperature does not provide better therapeutic results!
- Leave the DC connector connected to the device at all times for **non-battery model**. Frequent connection and disconnection of the connector may result in the damage of the wire inside the insulation. The connector is connected deep into the wood for avoiding easy disconnection that may interrupt the session. If the connector has to be reconnected, make sure the connector gets inserted as deep as it is designed.
- For **battery model**, use only the dedicated charger included in the product. Using other chargers may result in battery damage or explosion. Full charging takes about 3 hrs and the operation time is between 4 to 7 hours depending on the temperature set up.